

LARISA FACE COVER S.A.
6th KLM NATIONAL ROAD LARISSA
ATHENS, 41110
GR

Client Account Number: A00909959EFO
Eurofins Quote Number: N3OPPT20008601

General Method Reference: EN 14683:2019/AC 2019 App D

Analysis	Specification	Result	Unit
----------	---------------	--------	------

Eurofins Sample Number N720AA0310-1			
3ply non woven disposable surgical masks - aliquot 1; Lot # L202891; Original Received Date 19-Oct-2020;			
Microbial Cleanliness (Bioburden)			
Aerobic count	----	0.80	CFU/g
Aerobic count	----	3.00	CFU/sample
Yeasts and Moulds	----	0.80	CFU/g
Yeasts and Moulds	----	3.00	CFU/sample
Total Bioburden	≤ 30	1.60	CFU/g
Total Bioburden	----	6.00	CFU/sample
Analysis Date: 19-Oct-2020 to 26-Oct-2020			

Eurofins Sample Number N720AA0310-2			
3ply non woven disposable surgical masks - aliquot 2; Lot # L202891; Original Received Date 19-Oct-2020;			
Microbial Cleanliness (Bioburden)			
Aerobic count	----	< 0.80	CFU/g
Aerobic count	----	< 3.00	CFU/sample
Yeasts and Moulds	----	2.39	CFU/g
Yeasts and Moulds	----	9.00	CFU/sample
Total Bioburden	≤ 30	3.19	CFU/g
Total Bioburden	----	12.00	CFU/sample
Analysis Date: 19-Oct-2020 to 26-Oct-2020			

Eurofins Sample Number N720AA0310-3			
3ply non woven disposable surgical masks - aliquot 3; Lot # L202891; Original Received Date 19-Oct-2020;			
Microbial Cleanliness (Bioburden)			
Aerobic count	----	1.51	CFU/g
Aerobic count	----	6.00	CFU/sample
Yeasts and Moulds	----	1.51	CFU/g
Yeasts and Moulds	----	6.00	CFU/sample
Total Bioburden	≤ 30	3.02	CFU/g
Total Bioburden	----	12.00	CFU/sample
Analysis Date: 19-Oct-2020 to 26-Oct-2020			

Eurofins Sample Number N720AA0310-4			
3ply non woven disposable surgical masks - aliquot 4; Lot # L202891; Original Received Date 19-Oct-2020;			

Analysis	Specification	Result	Unit
Microbial Cleanliness (Bioburden)			
Aerobic count	----	1.53	CFU/g
Aerobic count	----	6.00	CFU/sample
Yeasts and Moulds	----	3.06	CFU/g
Yeasts and Moulds	----	12.00	CFU/sample
Total Bioburden	≤ 30	4.59	CFU/g
Total Bioburden	----	18.00	CFU/sample
Analysis Date: 19-Oct-2020 to 26-Oct-2020			

Eurofins Sample Number N720AA0310-5			
3ply non woven disposable surgical masks - aliquot 5; Lot # L202891; Original Received Date 19-Oct-2020;			
Microbial Cleanliness (Bioburden)			
Aerobic count	----	1.60	CFU/g
Aerobic count	----	6.00	CFU/sample
Yeasts and Moulds	----	0.80	CFU/g
Yeasts and Moulds	----	3.00	CFU/sample
Total Bioburden	≤ 30	2.40	CFU/g
Total Bioburden	----	9.00	CFU/sample
Analysis Date: 19-Oct-2020 to 26-Oct-2020			

Sample Compliance Assessment
Samples N720AA0310-1 through N720AA0310-5 meet the requirement(s) for all listed test(s) where specifications were applied.

Contracted Company: Eurofins Cosmetic & Personal Care Italy (Vimodrone)
via Bruno Buozzi, 2, Vimodrone (MI), IT-20055 Italy InfoCosme@eurofins.com

Questions about this report should be directed to your project manager or the general email listed above.